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CELIMED® C-715
(model NB-216C)

CZ

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REV.: CELIMED C-715 (NB-216C)-03012025

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REV.: CELIMED C-715 (NB-216C)-03012025

TABLE OF CONTENT

INTENDED USE	63
SYMBOLS	63
SAFETY INFORMATION	64
PARTS AND COMPONENTS	67
SAFETY INFORMATION.....	68
USE THE DEVICE WITH AC POWER ADAPTER	68
PREPARING THE NEBULIZER FOR USE.....	69
USING THE DEVICE	69
CLEANING.....	71
CARING FOR THE DEVICE	72
STORING THE DEVICE	73
TROUBLESHOOTING	74
WARRANTY	75
CLASSIFICATION	75
PERIODIC SAFETY CHECKS.....	75
SPECIFICATIONS.....	75
IMPORTANT INFORMATION REGARDING ELECTROMAGNETIC COMPATIBILITY (EMC)	76
COMPLIANCE INFORMATION FOR EACH EMC TEST	77

INTENDED USE

Thank you for purchasing the CELIMED C-715 (NB-216C) compressor nebulizer.

Intended purpose: This compressor nebulizer is a medical device that used to aerosolize medications for inhalation by the patient to achieve the expected therapeutic effect. The device can be used with adult or pediatric patients.

Intended user: Healthcare professionals; patients or their caregivers.

Intended patient: Pediatric (defined by the prescribed medication) and adult patients

Applicable place: hospitals, clinics, home.

Contraindications: None.

EN

SYMBOLS

Symbols	Meaning
	Manufacturer
	EU-representative
	DISPOSAL PROCEDURE (Directive 2012/19 / EU-WEEE) Symbol for the marking of electrical and electronics devices according to Directive 2002/96/EC. The device, accessories and the packaging have to be disposed of waste correctly at the end of the usage. Please follow Local Ordinances or Regulations for disposal
 CLASS II	CLASS II
CE 0123	CE marking in conformity with EC directive 93/42/EEC
IP 22	Protected against solid foreign objects 12.5 mm in diameter and larger, and against drops of water when the housing is angled up to 15°
	Observe the instructions for use
 TYPE BF	Type BF Applied Part
	Medical device

SAFETY INFORMATION

Read all instructions carefully before use!

WARNING

Use the equipment only for its intended use as described in this manual.

- Use the equipment with medications only under the instruction of your physician
- The device, including accessories, should not be shared by two different persons.
- Do not use the equipment if it has any damaged parts, or it has fallen into water. Return the product to a service center for examination and repair.
- When there is a HF magnetic field, high-power equipment, air compressor, heating equipment or other equipment that may affect the use of the product, please confirm whether the nebulizer can work /fog out normally before use.
- Do not place the device in an environment that contains a lot of dust and fluff. Fluff in the air, large particles of dust or insects may block the filter port and prevent the device from fogging out properly. Therefore, check whether the filter port is blocked by foreign matter before use.
- When pets or children play, the air tube may be bent or the device may be accidentally dropped, affecting the product performance. Therefore, please protect the device to prevent pets and children from affecting the use of the product when use.
- The air tube is long, please properly place or store it in the storage bin of the nebulizer to prevent strangulation or suffocation risk to children.
- For regime of medication shall follow the instructions of your physician or licensed healthcare practitioner.
- Do not cover the compressor with a blanket, towel, or any other type of cover during using. This could result in the compressor overheating or malfunctioning.
- Do not use the device where the device may be exposed to flammable gas or vapors.
- Do not use only water in the nebulizer for nebulizing purposes.
- Always dispose of any remaining medication in the medication tank after each use. Use fresh medication each time you use the device.
- Do not leave the device or its parts where it will be exposed to extreme temperatures or changes in humidity, such as leaving the device in a vehicle during warm or hot months, or where it will be exposed to direct sunlight.
- Do not use the accessories from other manufacturers. Using such accessories may could have an impact on the device's safety and performance.

CAUTION

- Provide close supervision when this device is used by, on, or near infants, children or compromised individuals.

- Do not insert any object into the compressor.
- Make sure that the air filter is clean. If the air filter has changed color or has not been used for 60 days, replace the filter.
- Make sure that the nebulizer kit is correctly assembled, the air filter is properly installed, and the air tube is correctly connected to the compressor and the nebulizer kit. Air may leak from the air tube during use if not securely connected.
- Inspect the compressor (main unit) and the nebulizer parts each time before using the device. Make sure no parts are damaged, the nozzle and air tube are not blocked and the compressor operates normally.
- Do not use the device if the air tube is bent.
- Do not block the air filter cover.
- Do not alter the baffle, the nozzle in the medication tank or any part of the nebulizer kit.
- Do not add more than 10ml of medication to the medication tank.

EN

CAUTION

- The material of the device and accessories in contact with human body has been tested for biological safety, among which both the skin irritation and sensitization meet the standard requirements of ISO10993 series.
- Do not operate the device at temperatures greater than .
- Do not tilt the nebulizer kit so the angle of the kit is greater than . Medication may flow into the mouth.
- Do not shake the nebulizer kit while using the device.
- Do not subject the compressor, or any of the components to strong shocks, such as dropping on the floor.
- This device is approved for human use only.
- Do not disassemble or attempt to repair the device or components.
- Use the device only for its intended use as described in the instruction manual. Do not use attachments not recommended by the manufacturer.
- Dispose of the device, components and optional accessories according to applicable local regulations. Unlawful disposal may cause environmental pollution.
- Make sure that the air tube is securely attached to the compressor (main unit) and nebulizing parts, and does not come loose. Twist the air tube slightly when inserting it into the connectors to avoid the tube disconnecting during use.
- Do not use the device in the MRI environment.
- The patient should choose the suitable size of mask based on their face size.

RISK OF ELECTRICAL SHOCK

- Do not use the compressor (main unit) and the AC adapter while they are wet.


WARNING

- (a) Every time before use, please check the adaptor, once damaged, never to use.
- (b) Please clean the plug of adapter plug at least once a year. Too much dust on plug may cause the fire
- Do not plug or unplug the power cord into the electrical outlet with wet hands.
 - Do not immerse the compressor (main unit) in water or other liquid.
 - Do not spill water or other liquids on the compressor .These parts are not waterproof. If liquid spills on these parts, please unplug the power cord and wipe off the liquid with gauze or other soft absorbent material immediately.
 - Do not use or store the device in humid locations or outdoors. Use the device within the operating temperature and humidity.
 - Do not overload power outlets. Plug the device into the appropriate voltage outlet.
 - Do not use extension cords. Plug the power plug directly into the electrical outlet.
 - Unplug the power plug from the electrical outlet after using the device. Never leave this product unattended when plugged in.
 - Unplug the power plug from the electrical outlet before cleaning the device.
 - Completely read all of the instructions included the optional accessories before using them.

MAINTENANCE AND STORAGE

- Keep the device out of the reach of unsupervised infants and children. The device may contain small parts that can be swallowed.
- Do not leave the cleaning solution in the nebulizer parts. Rinse the nebulizer parts with clean hot tap water after disinfecting.
- Wash the nebulizer parts after each use. Dry the parts immediately after washing.
- Do not store the air tube with moisture or medication remaining in the air tube. This could result in infection as a result of bacteria.
- Store the device and the components in a clean, safe location.
- Do not carry or leave the nebulizer with medication in the medication tank.
- Do not place or attempt to dry the device, components or any of the nebulizer parts in a microwave oven.

PARTS AND COMPONENTS

1. Main unit
2. Screen
3. On/off button
4. Countdown button
5. Pressure adjust button
6. Rubber foot
7. Air plug
8. Dust Screen cover
9. Type-C port

10. Holder for Nebulizer Kit
11. Nebulizer Kit
12. Adult mask
13. Child mask
14. Mouthpiece
15. Air Tube
16. Air Filter
17. Adapter

SAFETY INFORMATION

The followings are maintenance and repair which can be taken by operator, or which must be operated by manufacturer or distributor.

service and maintenance	Responsible
change the inhalation tube	operator
change the applied part	operator
change the air filter	operator
clean the surface of the device	operator
daily cleaning and disinfecting	operator
All components which need to be repaired or changed by separating the device	Authorized service center of distributor or manufacturer

WARNING

- Do not modify this equipment without authorization of the manufacturer.
- Do not service or maintain the device while in use with the patient.

USE THE DEVICE WITH AC POWER ADAPTER

Use the AC adapter as the power supply

- Insert the AC adapter cord into the jack on the side of the device. Insert the AC adapter plug into the outlet. To remove the AC adapter, disconnect the adapter plug from the AC outlet first and then disconnect the cord from the device's jack.

Caution:

To avoid electric shock, use only the exclusive AC adapter that can be purchased from authorized dealers. Other adapter may damage the device.

When the AC adapter is used as an isolating means, the AC adapter plug shall be inserted into the outlet nearby the operator, making it easy to disconnect the device from the outlet.

If long time work, remove the plug after the adapter cools, and prevent burns.

Caution: Optional AC adapter technical feature:

Output voltage: 12 V

Max. output current: 1.35 A

PREPARING THE NEBULIZER FOR USE

CAUTION

- Clean and disinfect set of spray chamber, mouthpiece and inhalation mask before using them for the first time after purchase.
- If the device has not been used for a long period of time please cleans and disinfects the set of spray chamber, mouthpiece and inhalation mask before using them.
- To prevent the risk of cross infection, different persons should never use the same nebulizer kit.



USING THE DEVICE

1. At normal use, place the device on a stable, sturdy and flat surface horizontally, such that the unit can be easily reached when you are seated.


EN

Caution: place the device at least 10 cm distance from walls.


2. Rotate the inhalation top counterclockwise to remove the inhalation top from the medication tank.
3. Add the correct amount of prescribed medication to the medication tank.
4. Turn the inhalation top clockwise until securely closed
5. Attach the desired inhalation accessory.
6. Connect the inhalation hose to the device and to its other end the chosen accessory (child's inhalation mask, adult mask or mouthpiece).
7. Remove the power source/USB adapter from the package and connect it to the side of the device/inhaler (see Fig. 9, on page 67 with a description of the device).
8. Insert the AC adapter plug into the outlet, the main unit powered on, entering standby mode, the screen will always display as shown in Figure 1. Press the on/off button  to start the device, „OFF“ will turn off, „ON“ will be on, and the device begins to run. The time display starts counting from 00 (as shown in Figure 2). During this period, press the on/off button  again, and the device will stop working and enter the standby mode (as shown in Figure 3). If there is no operation after working for 30 minutes, the „30“ will flash three times and then the device will stop running. The display will return to standby mode (as shown in Figure 3).

9. Pressure adjust function: In standby mode of the device, the initial pressure of the device is the small gear, and pressing the pressure adjust button  can adjust the pressure. There are three gears in total (gear changes are shown in Figure 4). The pressure level can also be adjusted during the device is running (as shown in Figure 5).


10. Countdown functional In standby mode, press the countdown button  to turn on the countdown function. „C/D“ remains on, „00“ flashes (as shown in Figure 6).

(1) “00” will flash 5 times and stop flashing if there is no operation . The device will return to standby mode. (as shown in Figure 7).

(2) Within „00“ flashing 5 times, press the countdown button to set the countdown time from 5 to 30 minutes. The smallest change unit is 5 minutes (as shown in Figure 8). Countdown time will flash 5 times before you finish the set. Countdown time is confirmed when the time number stops flashing, and the device starts running. (as shown in Figure 9). The display will start countdown.

During the operation of the device, you can press the countdown button to stop the device and re-enter the countdown setting (Figure 10).

USING THE CHILD MASK or THE ADULT MASK

Place the mask over the nose and mouth. Pull the elastic strap over the head. Gently pull on the strap to securely hold the mask over the nose and mouth. Inhale the medication. Exhale normally through the mask.

CLEANING

WARNING

- Ensure that the remove the plug before the clean and disinfection.
- Do not place or attempt to dry the device, components or any of the nebulizer parts in a microwave oven.

DURABLE PERIOD

As with any plastic parts, nebulizer kit and their accessories are affected by a certain amount of wear and tear when used and hygienically prepared on a frequent basis. Over time, this can lead to a change in the aerosol, which can have a negative effect on the efficiency of the treatment.

Durable period:

Compressor (main unit) 5 years

Nebulizer set and other accessories 300 inhalation cycles/use

Air filter 60 days

EN

PREPARATION BEFORE CLEANING

The nebulizer kit and the accessories used must be cleaned thoroughly after each application and disinfected at least once a day.

1. Remove the accessory (mask or mouthpiece) from the nebulizer kit.
2. Detach the air tube from the nebulizer kit.
3. Dismantle the nebulizer kit into its individual parts.
4. Make sure that all medication residues are removed from the nebulizer kit.
5. Rinse all parts of the nebulizer cup and accessory under running tap water. This must be done especially thoroughly if cleaning and disinfection are not carried out straight away. The air tube does not need to be cleaned or disinfected. If condensation occurs in the air tube, please clear it as follows.

The air tube does not need to be cleaned or disinfected. If condensation occurs in the air tube, please clear it as follows.

1. Connect the air tube to the compressor.
2. Switch the compressor on.
3. Leave the compressor running until any condensation in the tubing has been removed by the air flowing into the tube.

CLEANING

1. Place all disassembled parts in warm tap water (40–50°C) for at least 10 min. If necessary, use a clean brush to remove loose dirt (the brush must be reserved exclusively for this purpose).
2. Rinse all parts thoroughly in running water (the water flow should more than 3 liters per minute) and each part for at least 1 min.

3. Whether disinfect immediately or not, dry all parts with a new medical gauze. You can remove excess water more quickly by shaking them.

CAUTION

- Inadequate disinfection encourages the growth of bacteria and thus increases the risk of infection. Adequate cleaning with disinfection can only be assured if application time are adhered to, and if all individual parts are completely immersed in the solution for the entire application time. There must not be any cavities or air bubbles.

DISINFECTING

After cleaning, disinfect the dismantled nebulizer kit and accessor yused (except for the air tube). Effective disinfection is only possible if the nebulizer kit and accessory have been cleaned. Use 2.0% (w/w) Hydrogen Peroxide solution to disinfect the nebulizer kit and accessories. To ensure safety when handling chemicals, follow the instructions for use of the disinfecting agent, particularly the accompanying safety instructions. For complete disinfection of nebulizer kit and accessory, thoroughly clean, rinse and rough dry surfaces before immersed in chemical disinfection solution.

1. Place all pre-cleaned, rinsed and dried disassembled components in undiluted 2.0% (w/w) Hydrogen Peroxide solution. Ensure that all components' surfaces have been fully contacted with solution.
2. Once all components have been immersed and all surfaces in contact with the disinfectant solution, soak all components for 8 minutes. Track the soak time using a timer.
3. Rinse all components with sterile water or potable tap water.
4. Dry the parts with a new medical gauze as soon as disinfection has finished. You can remove excess water more quickly by shaking them.
5. Dispose of the used solution.

CAUTION

- A damp environment may encourage the growth of bacteria. Therefore, dry the parts with a new medical gauze as soon as disinfection has finished. The risk of infection is reduced when the parts are dried.

CARING FOR THE DEVICE

To keep your device in the best condition and protect the unit from damage follow these directions: :

CLEANING THE COMPRESSOR

Clean the casing of the main unit by using a soft cloth moistened with water or a mild detergent. Do not use abrasive cleaners. Dry the casing immediately using a soft clean cloth.

CHANGE THE AIR FILTER

Change the air filter every 60 days even if the air filter does not appear dirty. If the air filter appears dirty, or if water or medication has spilled on the air filter, replace with a new air filter immediately.

1. Pull the air filter cover to remove from the front side of the compressor.
2. Remove the dirty filter with hand

CAUTION:

Do not attempt to wash or clean the air filter. Damp air filters can cause blockages. Do not substitute cotton or any other material for the air filter.

CAUTION:

Wash the air filter cover regularly to prevent any blockage in the cover. Do not boil. Make sure the cover is dry before inserting the new air filter.

3. Insert a new air filter into the air filter cover.

EN

CAUTION:

Before inserting the new air filter makes sure the air filter is clean and free of dust. Do not operate the device without the air filter.

4. Put the air filter cover back on the compressor.
5. Do not maintain or service the device while it is in using.

STORING THE DEVICE

1. Put the nebulizer kit and the inhalation accessory (mouthpiece, or mask) in a dry, sealed bag. Place the bag in the storage pocket of the carrying bag.
2. Coil the air tubing and place in the storage pocket.
3. Place the compressor in the carrying bag.
4. Zipper the bag closed.
5. Store the device in a safe, clean location. Do not store the device in extreme hot or cold temperature, high humidity or in direct sunlight.

WARNING:

- Do not leave the device or its parts where it will be exposed to extreme temperatures or changes in humidity, such as leaving the device in a vehicle during warm or hot months, or where it will be exposed to direct sunlight.
- Keep the device out of the reach of unsupervised infants and children. The device may contain small parts that can be swallowed.

CAUTION:

Do not carry or leave the nebulizer with medication in the medication tank. Do not disassemble or attempt to repair the device or components.

TROUBLESHOOTING

TROUBLESHOOTING GUIDE

PROBLEM	CAUSE	SOLUTION
No power on unit when the power switch is on.	The AC power cord is not plugged into an electrical outlet.	Turn the power switch off. Plug the power plug into an electrical outlet. Turn the device on.
No nebulization or low nebulization rate when the power is on.	No medication in the medication tank. Too much or too little medication in the medication tank.	Add the correct amount of prescribed medication to the medication tank.
	The nebulizer kit is not correctly assembled.	Make sure the nebulizer kit is correctly assembled and the inhalation accessory is correctly attached.
	The nebulizer kit is tilted at an incorrect angle.	Hold the nebulizer kit correctly. Do not tilt the nebulizer kit so the angle of the kit is greater than 45 degrees.
	The air tube is incorrectly attached.	Make sure the air tube is correctly attached to the compressor and the nebulizer kit.
	The air tube is folded or damaged. The air tube is blocked.	Make sure the air tube is not folded, kinked or bent. Inspect the air tube for any damage. Replace the air tube if damaged.

The device is very hot.	The compressor is covered. The device has been used for longer than 20 minutes.	Do not cover the compressor with any type of cover during use. Turn the device off. Wait 40 minutes before using the device again.
-------------------------	---	--

WARRANTY

The warranty for the main unit is two years. The warranty does not cover the consumer material (accessories and consumer parts). Warranty will not be accepted if tempered with the device by using it in any other way than described by the instruction manual, or if a non-authorized service tries to open and / or repair the device.

EN

CLASSIFICATION

Equipment Classification with respect to protection from electric shock: CLASS II

Degree of protection from electric shock: TYPE BF

Degree of protection against ingress of water is rated as IPX2.

Equipment not suitable for use in the presence of flammable anesthetic mixture with air, oxygen or nitrous oxide.

PERIODIC SAFETY CHECKS

If you use the device with power adapter, preventive inspection and maintenance to be performed including the frequency of such maintenance.

- Every time before use, please check the adapter, once damaged, never to use.
- Please clean the plug of adapter plug at least once a year. Too much dust on plug may cause fire

SPECIFICATIONS

Product Category	Aerosol Therapy Equipment
Product Description	Compressor nebulizer
Model (code)	CELIMED C-715 (NB-216C)
Dimensions	166 × 136 × 104 mm
USB adapter/ IN	100–240 V AC, 50/60 Hz, 700 mA Max
USB adapter/ OUT	12V DC, 1,35 A

Extreme Pressure	30–60 psi/207–414 kPa
Free Flow Range	≥ 4 L/min.
Nebulizing Pressure	6 psi/42 kPa
Noise Level	< 55 dB
Max capacity of nebulizer kit	max. 10 ml
Particle Size	1~5 um ≥60%
Nebulizing Rate	0.3–0.5 ml/min
Operating Conditions	+5 °C–+40 °C, 15%–90% RH, 700~1060 hPa
Transport & Storage Conditions	-25 °C–+70 °C, 15%–90% RH, 700~1060 hPa
Protection against injury caused by electric shock	Category II
High Altitudes (m)	2000 m
Contents:	Compressor, USB power adapter+cable, spray chamber set, inhalation hose (PVC), mouthpiece, mask for adults (PVC), children's mask (PVC), spare air filters, instruction manual.

NOTES:

- Subject to technical modification without prior notice.
- Please note that specifications may vary with medication type used.
- Do not use the device where it may be exposed to flammable gas.
- This unit confirms to EMC Standard IEC60601-1. However, if it is used together with other medical devices or electrical equipment, they may influence the operations of one of devices. Please follow any instructions in the manuals and use all devices correctly.

IMPORTANT INFORMATION REGARDING ELECTROMAGNETIC COMPATIBILITY (EMC)

Warning:

- Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the equipment, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

- Under the test conditions specified in immunity, the product can provide the basic safety and essential performance.
- The Emissions characteristics of this equipment make it suitable for use in home healthcare environment (CISPR 11 Class B).

COMPLIANCE INFORMATION FOR EACH EMC TEST

Electromagnetic Emission (Home Healthcare Environment)	
Emission test (IEC 60601-1-2:2014)	Compliance
Conducted and radiated RF emissions	CISPR 11 Group 1 Class B
Harmonic emissions IEC 61000-3-2	Class A
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies

EN

Declaration-Electromagnetic Immunity (Home Healthcare Environment)		
Immunity test	EC 60601 test level	Compliance level
Conducted RF IEC 61000-4-6:2013	3V 150 kHz to 80 MHz 6 V in ISM and between 0.15 MHz and 80 MHz	3V 150 kHz to 80 MHz 6 V in ISM and between 0.15 MHz and 80 MHz
Radiated RF IEC 61000-4-3: 2006+A1: 2007+A2: 2010	10 V/m 80 MHz až 2,7 GHz	10 V/m
Electrostatic discharge (ESD) IEC 61000-4-2: 2008	±8 kV contact±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact±2 kV, ±4 kV, ±8 kV, ±15 kV air
Electrical fast transient/ burstIEC 61000-4-4: 2012	±2 kV for power supply lines	±2 kV for power supply lines
Surge IEC 61000-4-5: 2005	±0.5 kV, ±1 kV line(s) to lines	±0.5 kV, ±1 kV line(s) to lines
Voltage dips, short interrup- tions and voltage variations on power supply input linesIEC 61000-4-11: 2004	0% UT, 0.5 Cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°0% UT, 1 Cycle and 70% UT, 25/30 cycles- single phase: at 0°0% UT, 250 cycles	0% UT, 0.5 Cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°0% UT, 1 Cycle and 70% UT, 25/30 cycles- single phase: at 0°0% UT, 250 cycles

NOTE: The EUT is the a.c. mains voltage prior to application of the test level. The following phenomenon still fulfills the requirement of basic safety and essential performance.

* UT: 230 V~/50 Hz, the pressure of the EUT is deviation the normal value but the value is still more than 10 psi when flow is 4.5 L/min.

** UT: 230 V~/50 Hz, the EUT stop working when adding 0% UT, but the EUT can restore its normal mode automatically.

Declaration-IMMUNITY to proximity fields from RF wireless communications equipment

Immunity test	Testfrequency	Modulation	Maximum power	Immunity level	Compliance level
Radiated RFIEC 61000-4-3: 2006+A1: 2007+A2: 2010	385 MHz	**Pulse Modulation: 18 Hz	1,8 W	27 V/m	27 V/m
	450 MHz	*FM+5 Hz deviation: 1 kHz sine	2 W	28 V/m	28 V/m
	710 MHz 745 MHz 780 MHz	**Pulse Modulation: 217 Hz	0,2 W	9 V/m	9 V/m
	810 MHz 870 MHz 930 MHz	**Pulse Modulation: 18 Hz	2 W	27 V/m	27 V/m
	1720 MHz 1845 MHz 1970 MHz	**Pulse Modulation: 217 Hz	2 W	28 V/m	28 V/m
	2450 MHz	**Pulse Modulation: 217 Hz	2 W	28 V/m	28 V/m
	5240 MHz	**Pulse Modulation: 217 Hz	0,2 W	9 V/m	9 V/m
	5500 MHz				
	5785 MHz				

Note *: As an alternative to FM modulation, 50% pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

Note **: The carrier shall be modulated using a 50% duty cycle square wave signal.

When any serious incident that have occurred in relation to the device please contact the manufacturer or EC representative. Or when any failures occurred regarding device, please contact the manufacturer. The manufacturer reserves the right to make technical changes without notice in the interest of progress. Prior notices will not be given in case of any amendments within this manual. The mentioned trademarks and names are owned by the corresponding companies.

<p>Manufacturer:</p> 	<p>HONSUN (NANTONG) CO., LTD. No.8, Tongxing Road, Economic&Technical Development Area, 226 009 Nantong City, Jiangsu, P. R. China</p>
<p>EU-representative:</p> 	<p>Shanghai International Holding Corp. GmbH (EUROPE) Eiffestrasse 80, 20537 Hamburg, GERMANY</p>
<p>Importer for the Czech Republic and authorized service:</p> 	<p>Celimed s.r.o. Sociální péče 3487/5a, 400 11 Ústí nad Labem, Česká republika Pozáruční servis – tel.: 475 208 180 e-mail: info@celimed.cz www.celimed.cz</p>