# LAICA<sup>®</sup>

# 2

ANNI DI **GARANZIA**  
YEARS **GUARANTEE**  
ANOS DE **GARANTIA**  
ANOS DE **GARANTIA**  
JAHRE **GARANTIE**  
ANS DE **GARANTIE**  
XPONIA **ΕΓΓΥΗΣΗ**  
ANI DE **GARANTIE**  
ROK **ZÁRUKA**  
ROK **ZÁRUKA**  
ÉV **GARANCIJA**  
LET **GARANCIJE**  
GODINA **GARANCIJE**

**LAICA S.p.A.**

Viale del Lavoro, 10  
36048 Barbarano Mossano (VI)  
Italy  
Tel. +39 0444.795314  
info@laica.com  
Made in China

**www.laica.it**

Data - Date

--	--	--	--	--	--	--	--	--	--

Timbro rivenditore  
Retailer's stamp  
Sello del revendedor  
Carimbo do revendedor  
Stempel des Händlers  
Cachet du revendeur  
Σφραγίδα καταστήματος  
Stampilă distribuitor  
Razítko prodejce  
Pečiatka predajcu  
Bolti pecsét  
Pečät prodajalca  
Pečät prodavaoca

**LAICA<sup>®</sup>**  
**TAGLIANDO**  
**DI CONTROLLO**

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Prodotto da: Vapo Healthcare Co., Ltd. Southern unit of third floor, building B, No 99 Yudai West Rd, High tech district, KunShan, Suzhou, Jiangsu Province, China

EC	REP	MedNet GmbH, Borkstrasse 10, 48163 Münster, Germany
----	-----	---

**Distribuito da: Laica S.p.A.**  
**Viale del lavoro, 10 – 36048 Barbarano Mossano (VI) – Italy**  
**Phone: +39.0444.795314 – [www.laica.it](http://www.laica.it)**  
 Made in China

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**ULTRASOUND AEROSOL THERAPY DEVICE - TYPE VP-U1 Instructions and warranty**

Dear Customer, Laica would like to thank you for the preference accorded to this product, designed according to the criteria of reliability and quality in order to meet your complete satisfaction.

**IMPORTANT PLEASE READ CAREFULLY BEFORE USE KEEP FOR FUTURE REFERENCE**

The instruction manual must be considered an integral part of the product and must be kept for its entire life. If the appliance is transferred to another owner, its documentation must also be transferred in its entirety. To ensure the safe and correct use of the product, the user must carefully read the instructions and warnings contained in the manual insofar as they provide important information regarding safety, user and maintenance instructions. In the event of loss of the instruction manual or the need to receive more information or clarification, please contact us at the address below:  
Laica S.p.A. - Viale del Lavoro, 10 - 36048 Barbarano Mossano (VI) – Italy - Tel. +39 0444.795314 - info@laica.com - www.laica.it

This aerosol is an effective tool for home treatment of disorders of the upper and lower respiratory tract. It is compact, quiet, easy to use, built in compliance with the current European regulations regarding design criteria for the safety of the electrical appliances for medical use ( Dir.93/42/EEC and Standard EN13544-1).

**KEY TO SIMBOLS**

Warning

Prohibition

Caution!  
Read the user instructions carefully

Symbol for "Type BF applied parts"

Complies with Medical Device Directive 93/42/EEC

Manufacturer

Production date

European Representation

Serial number

Batch production number

Class II appliance

**IP22:**

Degree of protection of electrical equipment covers, where the first figure indicates the degree of protection against penetration by solid foreign bodies (from 0 to 6) and the second figure the degree of protection against penetration by liquids (from 0 to 8).

**SAFETY WARNINGS**

Prior to use, ensure the product is intact and that there is no visible damage. If in doubt, contact the retailer.

- Keep the plastic bag away from children: danger of suffocation.
- This device must be used only for its intended purpose and according to the user instructions. All other types of use shall be considered inappropriate and therefore hazardous. The manufacturer shall not be deemed liable for any damage caused by inappropriate or incorrect use.
- People aged 14 years or older, people with physical, sensory or mental disabilities or unskilled people, may use and maintain the product only under adult supervision. Children must not play with the device.
- Ensure that the operator, i.e. the patient, is able to use the appliance correctly.
- **Keep all parts of the packaging as well as the air hose and adaptor away from children (risk of suffocation and strangling)**
- Before connecting the device, make sure the mains voltage corresponds to the voltage on the plate affixed to the product.
- ⊘ DO NOT install the appliance in areas that are difficult to access, for example too close to the wall as it would be too difficult to disconnect the device from the mains.
- ⊘ NEVER leave the device operating unsupervised. Once finished using, turn it off and unplug it from the mains.
- Handle the product with care, protect it against accidental impact, extreme temperature fluctuations, humidity, dust, direct sunlight and sources of heat.
- In the event of a failure and/or malfunction, turn off the device without tampering with it. Always contact the retailer for repairs.
- ⊘ DO NOT service the device during use.
- Be sure to have dry hands before connecting or disconnecting the device and when pressing the "U" switch.

- ⊗ DO NOT pull the power cord to unplug the adapter from the power outlet.
- ⊗ DO NOT use extensions, unwind the power cord and keep it away from sources of heat.
- ⊗ DO NOT tamper with the power cord for any reason. If damaged, contact your retailer.
  - Always use the "⏻" switch to turn the device off and then unplug the adapter.
  - Unplug the power cord and the adapter immediately after use and before cleaning.
- ⊗ DO NOT immerse the device and DO NOT pour water or other liquids onto it.
  - If the device falls in the water, do not try to pick it up. Immediately unplug it. DO NOT reuse the device after removing it from the water.
  - Only use original accessories of the manufacturer.
  - Before each use, make sure that the parts have been assembled correctly.
- ⊗ DO NOT use the device if the air tube is bent.
  - Make sure that the air filter is inserted correctly and that it is clean.
- ⊗ Do not use the device close to high-frequency electromagnetic emitters.
  - Keep the aeration opening free from dust and impurities. Avoid leaning the device on soft surfaces (i.e. bed or pillows).
  - The device is built with materials that could cause allergic reactions (air tube made of EVA, PVC masks, PP mouthpiece and PET medical ampoules).
- ⊗ DO NOT dry the device or its accessories in the microwave.

#### IMPORTANT NOTES FOR BEST OPERATION

- 1) Do not use the equipment if the nebuliser chamber is without water or if the water level is below the indicated level.
- 2) The disposal medical cup is **SINGLE USE**; replace it after each use to guarantee best nebulisation and to avoid damage to the equipment.
- 3) Do not use oily medicines with a viscosity over 3 cPs (centipoise), the unit to measure viscosity, for example cortisone, as it could damage the transducer.
- 4) Do not use essential oils (such as eucalyptus essence), active ingredients with solid elements (suspended) and viscous solutions as they could damage the transducer.
- 5) **After every use**, the nebuliser chamber (no.8 exploded view) and the transducer (no.9 exploded view) should be cleaned with a cloth dipped in pure alcohol.
- 6) Wait at least 15 minutes from one use to the next.
- 7) Clean the top cover (no.5 exploded view) of the device handling the flow duct with care (no.6 exploded view).

#### PRODUCT DESCRIPTION (see fig.1)

- 1) Main unit
- 2) "⏻" key

- 3) Operation warning light ("ON/OFF")
- 4) Water tank
- 5) Upper cover
- 6) Stream conveyor
- 7) Disposable medicine cup (8ml) with spare parts
- 8) Nebulization chamber
- 9) Transducer
- 10) Air connection pipe
- 11) Masks connector
- 12) Mask for adults
- 13) Mask for children
- 14) Mouthpiece
- 15) Adaptor plug
- 16) Adaptor

#### TECHNICAL SPECIFICATIONS

- Product name: Ultrasonic aerosol
- Trade name: MD6026P
- Applied part: masks and mouthpiece, BF type
- Installation classification and use:
  - main unit: portable,
  - attached part: hand-held,
  - No AP or APG
- Power supply: INPUT: 100-240V ~ 50/60Hz, 0.35A; OUTPUT: 12V, 1A
- Temporary use: 13 min. ON / 15 min. OFF.
- Life expectancy of the product in use (used to spray the drug twice a day for approximately 10 minutes): around 2 years for the main unit, about 180 hours for the accessories. The frequent use of the device may shorten the life expectancy of the product.
- Ultrasonic frequency: 1.7 MHz
- Spraying capacity: approximately 0.4 ml/minute
- MMAD: ≤ 5 micron
- Ampoule capacity: max filling volume 8 ml

#### Conditions during operation:

- Temperatures: min. 5°C - max. 40°C, with air humidity from 15% to 90% in the absence of condensation and atmospheric pressure from 700 to 1060hPa

#### Storage conditions:

- Temperature: min. -20°C - max. 70°C, with relative air humidity from 0 to 90% in the absence of condensation and with atmospheric pressure from 700 to 1060hPa

#### INSTRUCTIONS FOR USE

- 1) Before each use, carefully clean and disinfect the water tank, the cover, the masks connectors, the air connection tube, the mouthpiece and the masks according to paragraph "Cleaning and maintenance".
- 2) Place the device on a flat surface.
- 3) Remove the cover. The tank holds 15 ml. Fill with tap water up to the red level indicator. This water is used as liquid conducting to the drug the ultrasonic waves and is never nebulized. Do not put the tank directly under the tap to fill it with water. It is strictly forbidden to use salt solution and water with any chemical element. **Do not pour drugs directly in the tank.**
- 4) Fill the nebulizer with the drug, then place it in the water tank according to figure 3. Pour the drug paying attention not to exceed the maximum capacity (8ml). Follow the instructions for use and the dosing of the medicinal product you are using.
- 5) Close the tank with the proper cover.
- 6) Now, insert the pipe and the mouthpiece or the connector and the mask for adults/children in the cover, according to figure 6.

#### In case of inhalation through mouthpiece:

- In case of troubles at low respiratory system, for example cough or bronchial infections, use the mouthpiece.
- Sit upright.
  - Connect the air connection tube in the cover and in the other end insert the mouthpiece.
  - Tighten the mouthpiece well with the lips.
  - Inhale slowly and deeply through the mouth and exhale through the nose. To increase the efficacy of the therapy, hold your breath for a while after having inhaled. During the treatment, remain quiet and relaxed. Breathe regularly and not too quickly. If you wish to interrupt the inhalation shortly, remove the mouthpiece from the mouth and then reinsert it.

#### In case of inhalation with mask:

- The mask for inhalations is suitable for the treatment of the nasal and pharyngeal area.
- Sit upright.
  - Connect the air connection tube in the cover and in the other end insert the fitting and the mask.
  - Keep the mask lightly pressed against the nose; the mask has to cover mouth and nose without squashing them.
  - Inhale and exhale slowly and deeply through the nose. To increase the efficacy of the therapy, hold your breath for a while after having inhaled.

During the treatment, remain quiet and relaxed. Breathe regularly and not too quickly. If you wish to interrupt the inhalation shortly, remove the mask for a while and then place it again on nose and mouth.

- 7) Connect the transformer to the device and then to the power supply plug paying attention that the voltage is right.
- 8) Switch on the device pressing the "⏻" key. Make sure that the warning light switches on (blue colour). During inhalation, the device must compulsorily remain in vertical position: this assures its correct operation, as well as the maximum efficacy of therapy. Immediately interrupt the treatment if you feel unease during the aerosol-therapy.
- 9) The device has an automatic stop function, which deactivates the device when the medicine is almost exhausted or after 13 minutes of inhalation. The device can be switched off at any time by pressing the "⏻" switch. **IMPORTANT:** As with most equipment for aerosol therapy, at the end of the inhalation treatment, a certain amount of the pharmaceutical will remain in the nebulization ampoule: this is perfectly normal. This amount of pharmaceutical, also known as residual volume, cannot be nebulized.
- 10) Let the device cool down for some minutes, turn the device upside down or use an absorbing cloth to remove the residual liquid and perform the cleaning procedures according to paragraph "Cleaning and maintenance".
- 11) Wait at least 15 minutes from one use to the next.

#### REPLACEMENT OF THE PHARMACEUTICAL AMPOULE

The disposal medicine cup is **SINGLE USE**; replace it after each use to guarantee best nebulisation and to avoid damage to the equipment. Within the pack you will find ampoules which you can use for replacement purposes. Other ampoules can be purchased from your retailer.

#### MICROBIC CONTAMINATION

In the presence of conditions carrying the risk of microbial infection and contamination, personal use of the accessories only is advised (always consult a doctor).

#### MAINTENANCE AND DISINFECTION

- Attention! Hot surface! Before cleaning the appliance, allow it to cool down as the transducer heats up during use.
- After removing the plug, remove the cover and empty any residues from the tank.
- After each use, clean the spray chamber (no. 8) and the transducer (no. 9) with a cloth

soaked in ethyl alcohol. Clean the top cover (no. 5) of the appliance with a cloth soaked in water or vinegar so that no residue remains. When cleaning the lid, handle the flow conveyor carefully (no. 6).

- Clean the outside of the appliance using a soft cloth dampened with water and mild soap (do not clean the main unit by placing it under running water and prevent any liquid from entering the appliance).
- In case of daily use, it is advisable to disinfect or sterilise the accessories every day to limit the risk of infections. In this case, use a commercially available alcohol-based medical disinfectant respecting the doses and limitations provided by the disinfectant manufacturer. Immerse all the accessories in the disinfectant solution for the specified period. Completely remove the residues of the disinfectant by thoroughly rinsing the accessories under running water and allow them to air dry.
- ⊘ DO NOT use autoclave sterilisers, EO (ethylene oxide) gas sterilisers or low temperature plasma sterilisers.
- To disinfect the appliance it is also possible to use 70% alcohol (or an alcohol-based disinfectant spray). Check that no disinfectant residues remain in the appliance.
- The appliance, the tank and all the accessories must be perfectly dry before using them.
- Store the appliance in a cool, dry place.
- It is advisable to check the performance of the appliance every 2 years or after a repair. Contact the Laica assistance service (activity excluded from the warranty).

**DISPOSAL PROCEDURE (Dir. 2012/19/Eu-WEEE)**

 The symbol on the bottom of the device indicates the separated collection of electric and electronic equipment. At the end of life of the device, do not dispose it as mixed solid municipal waste, but dispose it referring to a specific collection centre located in your area or returning it to the distributor, when buying a new device of the same type to be used with the same functions. If the appliance to be disposed of is less than 25 cm, it can be returned to a retail location that is over 400 m² without having to purchase a new, similar device. This procedure of separated collection of electric and electronic devices is carried out forecasting a European environmental policy aiming at safeguarding, protecting and improving environment quality, as well as avoiding potential effects on human health due to the presence of hazardous substances in such equipment or to an improper use of the same or of parts of the same. Caution! The wrong disposal of electric and electronic equipment may involve sanctions.

**WARRANTY**

This equipment is covered by a 2-year warranty starting from the purchase date, which has

to be proven by means **of the stamp or signature of the dealer and by the fiscal payment receipt that shall be kept enclosed to this document.** This warranty period complies with the laws in force and applies only in case the consumer is an individual. Laica products are designed for domestic use only; therefore, their use for commercial purposes is forbidden. The warranty covers only production defects and does not apply, if the damage is caused by accidental events, misuse, negligence or improper use of the product. Use the supplied accessories only; the use of other accessories invalidate the warranty cover. Do not open the equipment for any reason whatsoever; the opening or the tampering of the equipment, definitely invalidate the warranty. Warranty does not apply to parts subject to tear and wear along the use, as well as to batteries, when they are included. Warranty decays after 2 years from the purchase date. In this case, technical service interventions will be carried out against payment. Information on technical service interventions, under warranty terms or against payment, may be requested by writing at info@laica.com. Repair or replacement of products covered by the warranty terms are free of charge. In case of faults, refer to the dealer; DO NOT send anything directly to LAICA. All interventions covered by warranty (including product replacement or the replacement of any of its parts) shall not be meant to extend the duration of the original warranty period of the replaced product.

The manufacturer cannot be made liable for any damages that may be caused, directly or indirectly, to people, things and pets resulting from the inobservance of the provisions contained in the suitable user manual, especially with reference to the instructions related to installation, use and maintenance of the equipment. Being constantly involved in improving its own products, Laica reserves the right to modify its product, totally or partially, with reference to any production need, without that this involves any responsibility whatsoever by Laica or its dealers.

 Manufactured by: Vapo Healthcare Co., Ltd. Southern unit of third floor, building B, No 99 Yudai West Rd, High tech district, KunShan, Suzhou, Jiangsu Province, China

 MedNet GmbH, Borkstrasse 10, 48163 Münster, Germany

**Distributed by: Laica S.p.A.**  
**Viale del lavoro, 10 – 36048 Barbarano Mossano (VI) – Italy**  
**Phone: +39.0444.795314 - www.laica.it**  
 Made in China

- DEVICE needs special precautions regarding electromagnetic compatibility (EMC) and needs to be installed and put into service according to the EMC information provided in the ACCOMPANYING DOCUMENTS;
- Portable and mobile RF communications equipment can affect DEVICE.
- Warning that the use of accessories, transducers and cables other than those specified with the exception of transducers and cables sold by the manufacturer of the DEVICE as replacement parts for internal components, may result in increased EMISSIONS or decreased IMMUNITY of the DEVICE.
- Warning that the DEVICE should not be used adjacent to or stacked with other equipment. During Voltage Dips test, the product will change its operation state from on to off, after test, it needs operator to press on/off button to recover. This phenomenon does not affect BASIC SAFETY or ESSENTIAL PERFORMANCE, according to the standard, it's deemed to comply with the requirement.

**A1 Electromagnetic Emissions-For all EQUIPMENT and SYSTEMS**

Guidance and manufacture's declaration electromagnetic emission		
The DEVICE is intended for use in the electromagnetic environment specified below. The customer or the user of the DEVICE should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The DEVICE uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emission CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	The DEVICE is suitable for use in all establishments, other than domestic and those directly connected to the public low-voltage power supply network that supplies building used for domestic purposes.
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

**A2 Electromagnetic Immunity -For Home Healthcare Environment EQUIPMENT and SYSTEMS.**

Guidance and manufacture's declaration electromagnetic immunity			
The DEVICE is intended for use in the electromagnetic environment specified below. The customer or the user of DEVICE should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	±8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floor are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines	±2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	30A/m	30A/m	
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	±1 kV line to line ±2 kV line to ground	±1 kV line to line ±2 kV line to ground	Mains power quality should be that of a typical commercial or hospital environment. If the user of the DEVICE requires continued operation during power mains interruptions, it is recommended that the DEVICE be powered from an uninterruptible power supply.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % UT; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°	0 % UT; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°	
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	The MANUFACTURER should consider reducing the minimum separation distance, based on RISK MANAGEMENT, and using higher IMMUNITY TEST LEVELS that are appropriate for the reduced minimum separation distance. Minimum separation distances for higher IMMUNITY TEST LEVELS shall be calculated using the following equation: $E = \frac{6}{d} \sqrt{P}$ Where P is the maximum power in W, d is the minimum separation distance in m, and E is the IMMUNITY TEST LEVEL in V/m.
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz	10 V/m	
NOTE UT is the a.c. mains voltage prior to application of the test level			

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**Distribuirna: Laica S.p.A.**